

A revolutionary new stereo microscope



6 Ergonomic design is not a 'nice to have' feature, but critical to operator comfort and productivity.





A revolutionary new stereo microscope

Since the invention of the microscope in the 1590's, there have been significant technological developments in the design of the optical microscope. However, one thing has remained constant in all this time:

Uncomfortable microscope eyepieces.

You may accept that it is just normal for microscope eyepieces to be a little uncomfortable, not very ergonomic and not that easy to use, **but it doesn't have to be this way ...**



A stereo microscope with comfortable eyepieces

The Ergo80 stereo microscope employs a revolutionary patented eyepiece design which does not require the user to precisely align their eyes with the eyepieces.

In fact, the user can sit back from the microscope (up to 38mm) and still view the image, meaning that uncomfortable and difficult-to-use eyepieces are a thing of the past.





Ergonomic design is not a 'nice to have' feature, but critical to operator comfort and productivity.







Patented eyepieces increases head freedom, reducing fatigue

Improving operator ergonomics is not just about improving comfort

Businesses choose Vision Engineering's ergonomic microscopes because they know their operators are more efficient, more accurate and more productive. So the operator benefits and so does the business.

Give your stereo microscope a health check!



An ergonomic body position makes the Ergo80 more comfortable, less fatiguing and, more importantly, much easier to use. Additionally, optimal operator ergonomics minimises the risk of repetitive strain-related injuries. A happy worker is a productive worker.

Freedom of head movement

An additional benefit of Vision Engineering's patented eyepieces is that users do not need to precisely align their eyes with the eyepieces. This freedom of movement reduces neck and back strain associated with the fixed body position of conventional microscope eyepieces.

A natural view of the subject

With conventional microscope eyepieces, operators must position their eyes very close to the eyepieces, blocking out ambient light. The intense light exiting the eyepieces causes the pupils to contract. Constant contraction and expansion of the pupils is the main cause of eye fatigue with microscopes.


With the patented eyepieces of Ergo80, users sit back from the eyepieces (up to 38mm), allowing ambient light into the eyes. Additionally, the light exiting the eyepieces is spread over a larger area, providing a more natural view of the subject.

Ability to wear glasses

With Ergo80, operators do not need to remove their glasses (or safety glasses) to use the microscope.



critical for re-work, repair, dissection and other manipulation tasks. Sitting back from the eyepieces provides users with much better peripheral vision, so they can co-ordinate hands in a natural manner.





System details

Incorporating over 50 years of proven optical experience in a high quality infinity corrected stereo zoom microscope, the Ergo80 delivers superb performance, with exceptional optics.

With an 8:1 zoom ratio, the Ergo80 has a standard magnification range of 8x - 64x (128x max.) allowing fast and accurate viewing of all subjects.



A compact, modular design allows accessories to be added to the configuration, without loss of clarity or contrast through the common main objective.

Objective lens	Zoom Magnification Range	Working Distance
0.7x	5.6x - 44.8x	130mm
1.0x	8x - 64x	78mm
2.0x	16x - 128x	34mm



Accessories



- Image captureRange of digital / video camera options.
- Imaging software
 Software solutions are available for image.

Software solutions are available for image mark-up, dimensioning and documentation purposes. View the online demo at **www.visioneng.com/dimensionone**







- Ergowedge
 Allows the eyepieces to be tilted -5° to -25°, maintaining optimum ergonomics for all users.
- Double iris diaphragm
 Adjustable, to increase depth of focus and sample contrast.

Technical details

Optical

- Patented 'Expanded-Pupil' eyepieces providing 10 times greater head freedom than conventional binocular microscope eyepieces.
- Standard magnification range 8x 64x
- 8:1 zoom magnification ratio

Illumination

- 360° multi-point LED ringlight provides long-life surface illumination
- LED substage illumination (bench stand only)

Dimensions


Bench Stand

- A = 560mm max B = 300mm

Boom Mount

D = 660mm max.

E = 400mm F = 280mm



Quality, service & support

More about Vision Engineering

Vision Engineering has been designing and manufacturing ergonomic microscopes for over 50 years.

With a philosophy of design innovation, Vision Engineering holds world patents for a number of optical techniques which significantly improve microscope ergonomics.



To date, over 300,000 'eyepieceless' and 'expanded image' microscopes have been installed for both industry and life science applications.

ISO 9001:2008

Vision Engineering Ltd is certified for the quality management system ISO 9001:2008.



FM 557119

Service & support

Vision Engineering has a network of international offices throughout Europe, Asia and North America, supported by fully trained distributor partners. Full user training, service, and support is available, ensuring the highest levels of customer support is maintained.

Find out more at www.visioneng.com »





as well as a complete line of optical and video non-contact measuring systems.

Jean Louis MENEGON

Représentant et Conseiller Région Sud Ouest

Tel: 06 76 08 96 83

Email: jlouismenegon@aol.com www.menegon-metrologie.com

Visiter le site web :

www.menegon-metrologie.com