



DIAWIRELESS

Precision by Simplicity



Both wireless solutions are adjusted to the DIATEST programme.

Easy handling and the usual reliability of all DIATEST products are guaranteed

Precision by Simplicity!

Jean-Louis MENEGON
Représentant & conseiller Sud Ouest

Tél: **06.76.08.96.83**
mail: jlouismenegon@aol.com
site web: www.menegon-metrologie.com



DIAWIRELESS

Precision by Simplicity



Radio transmission

- ✓ Operating time up to 4,000 hours
- ✓ Up to 2,000,000 data transmissions with 2 coin cells
- ✓ Operating distance approx. 20 m in production, up to 200 m in free field
- ✓ Electronic probes: resolution 0.1 μm
- ✓ Extremely high precision: ca. 0.2 μm with linearization
- ✓ Electronic probes: travel 4 mm
- ✓ Measured value can be seen on the holder
- ✓ Configuration: nominal size mm/inch, auto-power-off
- ✓ Type of measurements: min, max, difference, average
- ✓ Free Software for data evaluation in Excel
- ✓ Easy to use



Bluetooth®

- ✓ Analogue and digital display of readings
- ✓ Up to 7 measurements can be displayed at the PC
- ✓ Up to 7 hours continuous data stream
- ✓ Resolution up to 0.1 μm
- ✓ Electronic probe: digital probe with 10 mm travel
- ✓ Type of measurements: static, max., min., difference
- ✓ User information: work piece good / not good, change of battery via coloured LED on holder
- ✓ Configuration: nominal size, tolerance limits, adjustable ratio, measurement direction, auto-power-off

DIATEST



Jean-Louis MENEGON
Représentant & conseiller Sud Ouest

Tél: **06.76.08.96.83**

mail: jlouismenegon@aol.com

site web: www.menegon-metrologie.com

